



1150 Empire Central Place, Suite 116; Dallas, TX 75247  
1-866-57-LASER (866-575-2737) / 1-877-53-LASER (877-535-2737)



## THE BTL 4000 SERIES

The line offers an effective solution to all your physiotherapy needs, featuring full range of parameters and a great variety of functions combined with easy-to-use operation. It comes in two lines.

BTL-4000 Professional  
BTL-4000 Optimal

### The summary of main features

- Large Choice of single-diode laser probes 30-400mW (red and infrared) and laser Cluster probes 200-1800mW that can be simultaneously connected
- Continuous and pulse laser therapy (frequency modulation 0 – 5000Hz)
- Duty factor 10 – 90%
- Acupuncture set for laser acupuncture
- Nogier and EAV frequencies
- Built-in battery
- Built-in therapeutic protocols and encyclopedia for: rehabilitation, orthopedics, sports medicine, dentistry, gynecology, dermatology, ENT, pediatrics and general practice.
- 50 customer designed therapeutic protocols
- Choice of optical attachments for ENT, gynecology, urology and dentistry.
- BTL-4110 Laser Professional - one channel therapeutic laser (two probes can be connected simultaneously)
- BTL-4120 Laser Professional - two channel therapeutic laser (four probes can be connected simultaneously)



**BTL-4000 Series Introduction**

The BTL-4000 series includes electrotherapy units, ultrasound therapy units, therapy lasers and combi units



**Electrotherapy BTL-4000 Puls Professional**

The BTL-4000 Puls Professional electrotherapy products come as either 1 or 2-channel electrotherapy units.



**Ultrasound therapy BTL-4710 Sono Professional**

The BTL-4710 Sono Professional multifrequency therapeutic ultrasound.



**Laser therapy BTL-4000 Laser Professional**

The BTL-4000 Laser Professional low level laser therapy units



**Combi units BTL-4000 Professional**

The BTL-4000 Combi Professional units include two physiotherapy treatments in one unit.



**Electrotherapy BTL-4000 Puls Optimal**

The BTL-4000 Puls Optimal electrotherapy line features 1 and 2-channel electrotherapy units differing in the number of waveform types.



**Ultrasound therapy BTL-4710 Sono Optimal**

The BTL-4710 Sono Optimal therapy ultrasound unit.



**Laser therapy BTL-4110 Laser Optimal**

The BTL-4110 Laser Optimal low level laser therapy unit



**Combi units BTL-4000 Optimal**

The BTL-4000 Combi Optimal units include two physiotherapy treatments in one unit.



1150 Empire Central Place, Suite 116  
Dallas, TX 75247

1-866-57-LASER (866-575-2737)  
[www.MicrolightLaser.com](http://www.MicrolightLaser.com)

1-877-53-LASER (877-535-2737)  
[www.SportLaser.com](http://www.SportLaser.com)